

Bovine CDK5/Cyclin-dependent-like kinase 5 ELISA Kit

Catalog No: #EK19365

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CDK5/Cyclin-dependent-like kinase 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cyclin-dependent-like kinase 5
Target Name	Cyclin-dependent-like kinase 5
Other Names	Cell division protein kinase 5, Proline-directed protein kinase 33 kDa subunit, PDPK, Serine/threonine-protein kinase PSSALRE, Tau protein kinase II catalytic subunit, TPKII catalytic subunit, CDKN5
Accession No.	Q02399

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.053ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CDK5

Note: This product is for in vitro research use only and is not intended for use in humans or animals.