

# Mouse Cdk5/Cyclin-dependent-like kinase 5 ELISA Kit

Catalog No: #EK19366



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cdk5/Cyclin-dependent-like kinase 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent-like kinase 5
Target Name	Cyclin-dependent-like kinase 5
Other Names	CR6 protein kinase, CRK6, Cell division protein kinase 5, Serine/threonine-protein kinase PSSALRE, Tau protein kinase II catalytic subunit, TPKII catalytic subunit, Cdkn5, Crk6
Accession No.	P49615

## Application Details

Assay Type: Sandwich Detect Range: 12.5-800ng/mL Sensitivity: 6.34ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cdk5

Note: This product is for in vitro research use only and is not intended for use in humans or animals.