

Rat C1qa/Complement C1q subcomponent subunit A ELISA Kit

Catalog No: #EK19368



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat C1qa/Complement C1q subcomponent subunit A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Complement C1q subcomponent subunit A
Target Name	Complement C1q subcomponent subunit A
Accession No.	P31720

Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/mL Sensitivity: 1.67ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C1qa

Note: This product is for in vitro research use only and is not intended for use in humans or animals.