

# Rat Cd38/ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 ELISA Kit

Catalog No: #EK19370

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cd38/ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1
Target Name	ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1
Other Names	2'-phospho-ADP-ribosyl cyclase,2'-phospho-ADP-ribosyl cyclase/2'-phospho-cyclic-ADP-ribose transferase,2'-phospho-cyclic-ADP-ribose transferase,ADP-ribosyl cyclase 1,ADPRC 1,CD38H,Cyclic ADP-ribose hydrolase 1,cADPr hydrolase 1,CD38
Accession No.	Q64244

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.081ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name: Cd38

Note: This product is for in vitro research use only and is not intended for use in humans or animals.