

# Dog CCK/Cholecystokinins ELISA Kit

Catalog No: #EK19383



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog CCK/Cholecystokinins ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Cholecystokinins
Target Name	Cholecystokinins
Other Names	CCK
Accession No.	Q9TS44

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 7.9pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CCK

Note: This product is for in vitro research use only and is not intended for use in humans or animals.