

Pig CCK/Cholecystokinin ELISA Kit

Catalog No: #EK19385



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pig CCK/Cholecystokinin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Cholecystokinin
Target Name	Cholecystokinin
Other Names	CCK
Accession No.	P01356

Application Details

Assay Type: Sandwich Detect Range:0.312-20ug/mL Sensitivity:0.156ug/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CCK

Note: This product is for in vitro research use only and is not intended for use in humans or animals.