

Mouse Brca1/Breast cancer type 1 susceptibility protein homolog ELISA Kit

Catalog No: #EK19391



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Brca1/Breast cancer type 1 susceptibility protein homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Breast cancer type 1 susceptibility protein homolog
Target Name	Breast cancer type 1 susceptibility protein homolog
Other Names	RING-type E3 ubiquitin transferase BRCA1
Accession No.	P48754

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Brca1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.