

Human CD6/T-cell differentiation antigen CD6 ELISA Kit

Catalog No: #EK19393



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CD6/T-cell differentiation antigen CD6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human T-cell differentiation antigen CD6
Target Name	T-cell differentiation antigen CD6
Other Names	T12,TP120,CD6
Accession No.	P30203

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.08ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CD6

Note: This product is for in vitro research use only and is not intended for use in humans or animals.