

# Bovine ABCG2/ATP-binding cassette sub-family G member 2 ELISA Kit

Catalog No: #EK19397



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine ABCG2/ATP-binding cassette sub-family G member 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine ATP-binding cassette sub-family G member 2
Target Name	ATP-binding cassette sub-family G member 2
Other Names	Urate exporter, CD338
Accession No.	Q4GZT4

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/ml Sensitivity: 0.081ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ABCG2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.