

Rat Ambp/Protein AMBP ELISA Kit

Catalog No: #EK19401



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Ambp/Protein AMBP ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Protein AMBP
Target Name	Protein AMBP
Other Names	Itil
Accession No.	Q64240

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.404ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ambp

Note: This product is for in vitro research use only and is not intended for use in humans or animals.