

Human AGTR2/Type-2 angiotensin II receptor ELISA Kit

Catalog No: #EK19402



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human AGTR2/Type-2 angiotensin II receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Type-2 angiotensin II receptor
Target Name	Type-2 angiotensin II receptor
Other Names	Angiotensin II type-2 receptor,AT2
Accession No.	P50052

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.78ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AGTR2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.