

Mouse Cfi/Complement factor I ELISA Kit

Catalog No: #EK19404



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Cfi/Complement factor I ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Complement factor I |
| Target Name | Complement factor I |
| Other Names | C3B/C4B inactivator, If |
| Accession No. | Q61129 |

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.158ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cfi

Note: This product is for in vitro research use only and is not intended for use in humans or animals.