

Rat Cfi/Complement factor I ELISA Kit

Catalog No: #EK19405



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Rat Cfi/Complement factor I ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Rt |
| Specificity | Natural and recombinant Rat Complement factor I |
| Target Name | Complement factor I |
| Other Names | C3B/C4B inactivator, If |
| Accession No. | Q9WUW3 |

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.045ng/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cfi

Note: This product is for in vitro research use only and is not intended for use in humans or animals.