

Mouse Amacr/Alpha-methylacyl-CoA racemase ELISA Kit

Catalog No: #EK19406



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Amacr/Alpha-methylacyl-CoA racemase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-methylacyl-CoA racemase
Target Name	Alpha-methylacyl-CoA racemase
Other Names	2-methylacyl-CoA racemase,Macr1
Accession No.	O09174

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.2pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Amacr

Note: This product is for in vitro research use only and is not intended for use in humans or animals.