

Mouse Cd3g/T-cell surface glycoprotein CD3 gamma chain ELISA Kit

Catalog No: #EK19409



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cd3g/T-cell surface glycoprotein CD3 gamma chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-cell surface glycoprotein CD3 gamma chain
Target Name	T-cell surface glycoprotein CD3 gamma chain
Other Names	T-cell receptor T3 gamma chain, CD3g
Accession No.	P11942

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.056ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd3g

Note: This product is for in vitro research use only and is not intended for use in humans or animals.