

# Rat Cox5a/Cytochrome c oxidase subunit 5A, mitochondrial ELISA Kit

Catalog No: #EK19412



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cox5a/Cytochrome c oxidase subunit 5A, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cytochrome c oxidase subunit 5A, mitochondrial
Target Name	Cytochrome c oxidase subunit 5A, mitochondrial
Other Names	Cytochrome c oxidase polypeptide Va
Accession No.	P11240

## Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.312ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cox5a

Note: This product is for in vitro research use only and is not intended for use in humans or animals.