

Mouse Csnk2a1/Casein kinase II subunit alpha ELISA Kit

Catalog No: #EK19415



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Csnk2a1/Casein kinase II subunit alpha ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Casein kinase II subunit alpha
Target Name	Casein kinase II subunit alpha
Other Names	CK II alpha, Ckiia
Accession No.	Q60737

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.089ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Csnk2a1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.