

# Rat Csnk2a1/Casein kinase II subunit alpha ELISA Kit

Catalog No: #EK19416



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Csnk2a1/Casein kinase II subunit alpha ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Casein kinase II subunit alpha
Target Name	Casein kinase II subunit alpha
Other Names	CK II alpha
Accession No.	P19139

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.21ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Csnk2a1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.