

Mouse A2m/Alpha-2-macroglobulin ELISA Kit

Catalog No: #EK19418



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse A2m/Alpha-2-macroglobulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-2-macroglobulin
Target Name	Alpha-2-macroglobulin
Other Names	Alpha-2-M,Pregnancy zone protein,Pzp
Accession No.	Q61838

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/ml Sensitivity:39pg/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:A2m

Note: This product is for in vitro research use only and is not intended for use in humans or animals.