

Chicken ALB/Serum albumin ELISA Kit

Catalog No: #EK19423



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chicken ALB/Serum albumin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Serum albumin
Target Name	Serum albumin
Other Names	Alpha-livetin
Accession No.	P19121

Application Details

Assay Type: Sandwich

Detect Range:15.6-1000 pg/mL Sensitivity:21pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul

Assay Time:3-5h Detection wavelength:450 nm

gene_full_name:ALB

Note: This product is for in vitro research use only and is not intended for use in humans or animals.