

# Human ABCA4/Retinal-specific ATP-binding cassette transporter ELISA Kit

Catalog No: #EK19425



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ABCA4/Retinal-specific ATP-binding cassette transporter ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Retinal-specific ATP-binding cassette transporter
Target Name	Retinal-specific ATP-binding cassette transporter
Other Names	ATP-binding cassette sub-family A member 4,RIM ABC transporter,RIM protein,Stargardt disease protein,ABCR
Accession No.	P78363

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.9pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ABCA4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.