

Human ABCC3/Canalicular multispecific organic anion transporter 2 ELISA Kit

Catalog No: #EK19426

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ABCC3/Canalicular multispecific organic anion transporter 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Canalicular multispecific organic anion transporter 2
Target Name	Canalicular multispecific organic anion transporter 2
Other Names	ATP-binding cassette sub-family C member 3, Multi-specific organic anion transporter D, MOAT-D, Multidrug resistance-associated protein 3, CMOAT2, MLP2, MRP3
Accession No.	O15438

Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/mL Sensitivity: 33 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100 µl Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ABCC3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.