

# Rat Abcc3/Canalicular multispecific organic anion transporter 2 ELISA Kit

Catalog No: #EK19427

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Abcc3/Canalicular multispecific organic anion transporter 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Canalicular multispecific organic anion transporter 2
Target Name	Canalicular multispecific organic anion transporter 2
Other Names	ATP-binding cassette sub-family C member 3,MRP-like protein 2,MLP-2,Multidrug resistance-associated protein 3,Cmoat2,Mlp2,Mrp3
Accession No.	O88563

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.062ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Abcc3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.