

Mouse Abcc6/Multidrug resistance-associated protein 6 ELISA Kit

Catalog No: #EK19428



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Abcc6/Multidrug resistance-associated protein 6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Multidrug resistance-associated protein 6
Target Name	Multidrug resistance-associated protein 6
Other Names	ATP-binding cassette sub-family C member 6,Mrp6
Accession No.	Q9R1S7

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.062ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Abcc6

Note: This product is for in vitro research use only and is not intended for use in humans or animals.