

# Human ACADM/Medium-chain specific acyl-CoA dehydrogenase, mitochondrial ELISA Kit

Catalog No: #EK19435



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ACADM/Medium-chain specific acyl-CoA dehydrogenase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Medium-chain specific acyl-CoA dehydrogenase, mitochondrial
Target Name	Medium-chain specific acyl-CoA dehydrogenase, mitochondrial
Other Names	MCAD
Accession No.	P11310

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ACADM

Note: This product is for in vitro research use only and is not intended for use in humans or animals.