

Bovine ASPN/Asporin ELISA Kit

Catalog No: #EK19437



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine ASPN/Asporin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Asporin
Target Name	Asporin
Accession No.	Q3ZBN5

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.16ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ASPN

Note: This product is for in vitro research use only and is not intended for use in humans or animals.