

Human ASPN/Asporin ELISA Kit

Catalog No: #EK19438



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ASPN/Asporin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Asporin
Target Name	Asporin
Other Names	Periodontal ligament-associated protein 1,PLAP-1,PLAP1,SLRR1C,UNQ215/PRO241
Accession No.	Q9BXN1

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.056ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ASPN

Note: This product is for in vitro research use only and is not intended for use in humans or animals.