

# Rat Cckbr/Gastrin/cholecystokinin type B receptor ELISA Kit

Catalog No: #EK19439



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cckbr/Gastrin/cholecystokinin type B receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Gastrin/cholecystokinin type B receptor
Target Name	Gastrin/cholecystokinin type B receptor
Other Names	CCK-B receptor,Cholecystokinin-2 receptor,CCK2-R
Accession No.	P30553

## Application Details

Assay Type: SandwichDetect Range:78-5000pg/ml Sensitivity:40.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Cckbr

Note: This product is for in vitro research use only and is not intended for use in humans or animals.