

Mouse Apaf1/Apoptotic protease-activating factor 1 ELISA Kit

Catalog No: #EK19440



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Apaf1/Apoptotic protease-activating factor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apoptotic protease-activating factor 1
Target Name	Apoptotic protease-activating factor 1
Other Names	APAF-1
Accession No.	O88879

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Apaf1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.