

Rat Aifm1/Apoptosis-inducing factor 1, mitochondrial ELISA Kit

Catalog No: #EK19443



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Aifm1/Apoptosis-inducing factor 1, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Apoptosis-inducing factor 1, mitochondrial
Target Name	Apoptosis-inducing factor 1, mitochondrial
Other Names	Programmed cell death protein 8,Aif,Pdcd8
Accession No.	Q9JM53

Application Details

Assay Type: Sandwich Detect Range:31.2-2000pg/mL Sensitivity:14.3pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aifm1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.