

Mouse Adrb1/Beta-1 adrenergic receptor ELISA Kit

Catalog No: #EK19444



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Adrb1/Beta-1 adrenergic receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta-1 adrenergic receptor
Target Name	Beta-1 adrenergic receptor
Other Names	Beta-1 adrenoreceptor, Beta-1 adrenoceptor, Adrb1r
Accession No.	P34971

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.17ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adrb1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.