

# Bovine ABCC1/Multidrug resistance-associated protein 1 ELISA Kit

Catalog No: #EK19445



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine ABCC1/Multidrug resistance-associated protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Multidrug resistance-associated protein 1
Target Name	Multidrug resistance-associated protein 1
Other Names	ATP-binding cassette sub-family C member 1,Leukotriene C(4) transporter,LTC4 transporter,MRP1
Accession No.	Q8HXQ5

## Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.32ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ABCC1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.