

Human A1BG/Alpha-1B-glycoprotein ELISA Kit

Catalog No: #EK19448



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human A1BG/Alpha-1B-glycoprotein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Alpha-1B-glycoprotein
Target Name	Alpha-1B-glycoprotein
Other Names	Alpha-1-B glycoprotein
Accession No.	P04217

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.083ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:A1BG

Note: This product is for in vitro research use only and is not intended for use in humans or animals.