

Mouse A1bg/Alpha-1B-glycoprotein ELISA Kit

Catalog No: #EK19449



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse A1bg/Alpha-1B-glycoprotein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-1B-glycoprotein
Target Name	Alpha-1B-glycoprotein
Other Names	Alpha-1-B glycoprotein
Accession No.	Q19LI2

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.085ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: A1bg

Note: This product is for in vitro research use only and is not intended for use in humans or animals.