

Rat A1bg/Alpha-1B-glycoprotein ELISA Kit

Catalog No: #EK19450



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat A1bg/Alpha-1B-glycoprotein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-1B-glycoprotein
Target Name	Alpha-1B-glycoprotein
Other Names	Alpha-1-B glycoprotein, C44, Liver regeneration-related protein 1
Accession No.	Q9EPH1

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 36pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: A1bg

Note: This product is for in vitro research use only and is not intended for use in humans or animals.