

# Rat Bgn/Biglycan ELISA Kit

Catalog No: #EK19452



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Bgn/Biglycan ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Biglycan
Target Name	Biglycan
Other Names	Bone/cartilage proteoglycan I,PG-S1
Accession No.	P47853

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.086ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bgn

Note: This product is for in vitro research use only and is not intended for use in humans or animals.