

# Rat Cela2a/Chymotrypsin-like elastase family member 2A ELISA Kit

Catalog No: #EK19454



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cela2a/Chymotrypsin-like elastase family member 2A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Chymotrypsin-like elastase family member 2A
Target Name	Chymotrypsin-like elastase family member 2A
Other Names	Elastase-2, Elastase-2A, Ela2, Ela2a
Accession No.	P00774

## Application Details

Assay Type: Sandwich Detect Range: 3.12-200U/L Sensitivity: 1.34U/L Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cela2a

Note: This product is for in vitro research use only and is not intended for use in humans or animals.