

Human ATP6AP2/Renin receptor ELISA Kit

Catalog No: #EK19458

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ATP6AP2/Renin receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Renin receptor
Target Name	Renin receptor
Other Names	ATPase H(+)-transporting lysosomal accessory protein 2,ATPase H(+)-transporting lysosomal-interacting protein 2,ER-localized type I transmembrane adaptor,Embryonic liver differentiation factor 10,N14F,Renin/prorenin receptor,Vacuolar ATP synthase membrane sector-associated protein M8-9,ATP6M8-9,ATP6IP2,CAPER,ELDF10,HT028,MSTP009,PSEC0072
Accession No.	O75787

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ATP6AP2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.