

Human ADAD1/Adenosine deaminase domain-containing protein 1 ELISA Kit

Catalog No: #EK19459



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ADAD1/Adenosine deaminase domain-containing protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Adenosine deaminase domain-containing protein 1
Target Name	Adenosine deaminase domain-containing protein 1
Other Names	Testis nuclear RNA-binding protein, TENR
Accession No.	Q96M93

Application Details

Assay Type: Sandwich Detect Range: 6.25-400pg/mL Sensitivity: 2.1pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ADAD1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.