

Mouse Actn3/Alpha-actinin-3 ELISA Kit

Catalog No: #EK19463



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Actn3/Alpha-actinin-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-actinin-3
Target Name	Alpha-actinin-3
Other Names	Alpha-actinin skeletal muscle isoform 3,F-actin cross-linking protein
Accession No.	O88990

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.39ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Actn3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.