

# Mouse Bcl6b/B-cell CLL/lymphoma 6 member B protein ELISA Kit

Catalog No: #EK19468



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mouse Bcl6b/B-cell CLL/lymphoma 6 member B protein ELISA Kit         |
| Brief Description  | ELISA Kits   |
| Applications       | ELISA  |
| Species Reactivity | Ms   |
| Specificity        | Natural and recombinant Mouse B-cell CLL/lymphoma 6 member B protein |
| Target Name        | B-cell CLL/lymphoma 6 member B protein                               |
| Other Names        | Bcl6-associated zinc finger protein,Bazf                             |
| Accession No.      | O88282   |

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.084ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bcl6b

Note: This product is for in vitro research use only and is not intended for use in humans or animals.