

Human BCL6/B-cell lymphoma 6 protein ELISA Kit

Catalog No: #EK19469

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BCL6/B-cell lymphoma 6 protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human B-cell lymphoma 6 protein
Target Name	B-cell lymphoma 6 protein
Other Names	BCL-6,B-cell lymphoma 5 protein,BCL-5,Protein LAZ-3,Zinc finger and BTB domain-containing protein 27,Zinc finger protein 51,BCL5,LAZ3,ZBTB27,ZNF51
Accession No.	P41182

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.084ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BCL6

Note: This product is for in vitro research use only and is not intended for use in humans or animals.