

# Human ALK/ALK tyrosine kinase receptor ELISA Kit

Catalog No: #EK19470



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ALK/ALK tyrosine kinase receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human ALK tyrosine kinase receptor
Target Name	ALK tyrosine kinase receptor
Other Names	Anaplastic lymphoma kinase,CD246
Accession No.	Q9UM73

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 41.2pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ALK

Note: This product is for in vitro research use only and is not intended for use in humans or animals.