

## Human CDK16/Cyclin-dependent kinase 16 ELISA Kit

Catalog No: #EK19472

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CDK16/Cyclin-dependent kinase 16 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cyclin-dependent kinase 16
Target Name	Cyclin-dependent kinase 16
Other Names	Cell division protein kinase 16,PCTAIRE-motif protein kinase 1,Serine/threonine-protein kinase PCTAIRE-1,PCTAIRE1,PCTK1
Accession No.	Q00536

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.085ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CDK16

Note: This product is for in vitro research use only and is not intended for use in humans or animals.