

## Mouse Cdk16/Cyclin-dependent kinase 16 ELISA Kit

Catalog No: #EK19473

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cdk16/Cyclin-dependent kinase 16 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent kinase 16
Target Name	Cyclin-dependent kinase 16
Other Names	CRK5, Cell division protein kinase 16, PCTAIRE-motif protein kinase 1, Serine/threonine-protein kinase PCTAIRE-1, Crk5, Pctaire1, Pctk1
Accession No.	Q04735

## Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.46ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cdk16

Note: This product is for in vitro research use only and is not intended for use in humans or animals.