

Rat Bmp10/Bone morphogenetic protein 10 ELISA Kit

Catalog No: #EK19475



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Bmp10/Bone morphogenetic protein 10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bone morphogenetic protein 10
Target Name	Bone morphogenetic protein 10
Other Names	BMP-10
Accession No.	Q4AEG6

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.062ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bmp10

Note: This product is for in vitro research use only and is not intended for use in humans or animals.