

# Bovine ARG1/Arginase-1 ELISA Kit

Catalog No: #EK19476



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine ARG1/Arginase-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Arginase-1
Target Name	Arginase-1
Other Names	Liver-type arginase, Type I arginase
Accession No.	Q2KJ64

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 41.7pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ARG1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.