

# Rat Arg1/Arginase-1 ELISA Kit

Catalog No: #EK19477



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Arg1/Arginase-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Arginase-1
Target Name	Arginase-1
Other Names	Liver-type arginase, Type I arginase
Accession No.	P07824

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 12pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Arg1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.