

Human CENPB/Major centromere autoantigen B ELISA Kit

Catalog No: #EK19478



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CENPB/Major centromere autoantigen B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Major centromere autoantigen B
Target Name	Major centromere autoantigen B
Other Names	Centromere protein B,CENP-B
Accession No.	P07199

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.12ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CENPB

Note: This product is for in vitro research use only and is not intended for use in humans or animals.