

Mouse Cd52/CAMPATH-1 antigen ELISA Kit

Catalog No: #EK19482



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cd52/CAMPATH-1 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CAMPATH-1 antigen
Target Name	CAMPATH-1 antigen
Other Names	Lymphocyte differentiation antigen B7, CD52, Cdw52, Mb7
Accession No.	Q64389

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/ml Sensitivity: 0.083ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd52

Note: This product is for in vitro research use only and is not intended for use in humans or animals.